



		SIRE	
		R&W	B&W *RC
DAM	R&W	100% R&W	50% R&W 50% *RC
	B&W *RC	50% R&W 50% *RC	25% R&W 50% *RC
	B&W	0% R&W 100% *RC	0% R&W 50% *RC

(CDN, AUG\*14)

## VARIANT RED \*VRC

Variant Red Gene or Dominant Red is independent of the traditional red gene. It acts dominant over B&W.



OF PROGENY FROM \*VRC SIRE


(CDN, AUG\*14)

## TOP GTPI RED & WHITES



with GTPI  
greater than  
**2700**

CDCB APR\*21



Use **Elevate<sup>®</sup>**  
to confirm coat  
colour.

Use  
**SemexWorks<sup>™</sup>**  
to select for  
Red & White.

Use  
**Optimate<sup>™</sup>**  
to maximize  
your odds.

Red & White Sires With High Genetic Value